


## **Serotonin Transporter-Linked Polymorphic Region (5HTTLPR) and rs25531 SNP (M<sub>spl</sub>, L<sub>A</sub> / L<sub>G</sub>)**

**Serotonin Transporter (5HTT, Locus Symbol SLC6A4)**, which maps to 17q11.1-17q12 (Ramamoorthy et al., 1993), contains a **43 bp insertion/deletion** (ins/del, **5HTTLPR**) polymorphism in the 5' regulatory region of the gene (Heils *et al.*, 1996). It should be noted that due to an error in sequencing this was originally thought to be a 44 bp deletion (the highly repetitive nature of this site makes this a very excusable error). The ins/del in the promoter appears to be associated with variations in transcriptional activity: the long variant (**L**) has approximately three times the basal activity of the short promoter (**S**) with the deletion (Lesch *et al.*, 1996), although this is not a universal finding (Willeit et al., 2001, Kaiser et al., 2002). The **S** variant has been reported to be dominant over the **L** variant (Heils et al., 1996), although at least one report suggests that the **L** may be dominant over the **S** (Williams et al, 2003). Several investigators have reported that the 5-HTTLPR polymorphism affects serotonergic functions *in vivo*. Individuals with the **L/L** genotype were found to have significantly higher maximal uptake of serotonin into platelets compared to those with **L/S** or **S/S** genotypes (Nobile et al., 1999, Greenberg et al., 1999).

## Serotonin Transporter, SLC6A4 17q11.1-q12 minus strand



A depiction of the organization of the serotonin transporter showing the 5HTTLPR region and the positions of several SNPs that will be used in other analyses.  
Adapted from Heils et al, 1996 and Lesch et al, 1996.

The assay we use for the 5HTTLPR (Anchordoquy et al, 2003) is a modification of the method of Lesch *et al*, (1996). The primer sequences are from Gelernter et al. (1999).

**Forward: 5'- 6FAM - ATG CCA GCA CCT AAC CCC TAA TGT - 3',**  
**Reverse: 5'- GGA CCG CAA GGT GGG CGG GA - 3'.**

These primers yield amplicons of 419 (L) or 376 (S) bp.

## 5HTTLPR PCR Master Mix for 20 $\mu$ L reactions

(18  $\mu$ L Master mix + 2  $\mu$ L DNA)

Component	1 Tube vol ( $\mu$ L)	100 Tubes vol ( $\mu$ L)	Concentration of component in:		
			Stock	Master Mix	PCR
Water	9.3	930			
DMSO	2.0	200	100%	10.9 %	10%
10x Buffer II	2.0	200	10 x	0.109 x	1 x
MgCl <sub>2</sub>	1.6	160	25 mM	2.18mM	2.00 mM
dNTP+deazaGTP	2.0	200	2 mM (ea)	218 $\mu$ M	200 $\mu$ M (ea)
Forward	0.65	65	12 $\mu$ M	425 $\mu$ M	380 $\mu$ M
Reverse	0.65	65	12	425 $\mu$ M	380 $\mu$ M
AmpliTaq Gold®	0.2	20	5 units/ $\mu$ L	1 unit	1.0 unit
Total volume ( $\mu$ L)	18.4	1840			

### Preparation of dNTPs + 7-deaza-2-deoxy GTP

Component	volume ( $\mu$ L)	Concentration (mM)	
		Stock	Final
dATP	40	100	2
dTTP	40	100	2
dCTP	40	100	2
dGTP	20	100	1
deazaGTP	200	10	1
Water	1660		

### 5HTTLPR PCR Setup

Mastermix	18 $\mu$ L
DNA	1-2 $\mu$ L (20 ng or less)
Water	0-1 $\mu$ L
Total volume	20 $\mu$ L


### 5HTTLPR Touchdown PCR Cycling

1x	95 °C	10 min			
1x	95 °C	30 sec	65 °C	30 sec	72 °C 90 sec
1x	95 °C	30 sec	64 °C	30 sec	72 °C 90 sec
1x	95 °C	30 sec	63 °C	30 sec	72 °C 90 sec
1x	95 °C	30 sec	62 °C	30 sec	72 °C 90 sec
1x	95 °C	30 sec	61 °C	30 sec	72 °C 90 sec
1x	95 °C	30 sec	90 °C	30 sec	72 °C 90 sec
1x	95 °C	30 sec	59 °C	30 sec	72 °C 90 sec
1x	95 °C	30 sec	58 °C	30 sec	72 °C 90 sec
1x	95 °C	30 sec	57 °C	30 sec	72 °C 90 sec
1x	95 °C	30 sec	56 °C	30 sec	72 °C 90 sec
30x	90 °C	30 sec	55 °C	30 sec	72 °C 90 sec
1x	72 °C	30 min			
	4 °C	hold			

## 5HTTLPR Electrophoresis

2 µL PCR product  
 20 µL Hi-Di formamide  
 0.5 µL Genescan 500 (or 2500) Rox

Samples are analyzed on an ABI PRISM® 3130xl Genetic Analyzer using standard company protocols without modification



5HTTLPR. The figure above is reproduced from a run from an ABI PRISM® 3100 Genetic Analyzer. The amplicons are labeled with size in base pairs. The sizes given are those calculated from the DNA sequence, but in actual runs the size of the amplicon sizes calculated by the software are usually 3-4 bp greater. The figure above shows both alleles which generally have approximately equal peak heights. The red peaks are size standards (Genescan ROX 500).

The table below lists the frequencies of the two alleles in approximately 1000 subjects taken from the National Youth Survey Family Study. These are typical for a largely Caucasian population such as this.

Amplicon Size	Allele	Frequency
---------------	--------	-----------

376  
419

short  
long

.43  
.57

Notes:

For consistent results with this primer set the use of 10% DMSO and 7-deaza-2-deoxy GTP (Roche Applied Science, Indianapolis, IN) is essential.

Use a very good grade of DMSO. We use Sigma's Hybra-Max® grade or that supplied with New England Biolab's Phusion™ buffers.

We use touchdown PCR (Don et al, 1992) routinely as a simple short cut. It cuts down on the need to optimize annealing conditions for multiple primer sets when you want to do several loci in the same thermocycler.

There are many papers on this polymorphism, which can lead to confusion at first. Depending on the primer sets used and the nomenclature the authors use, the sizes of the reported long and short alleles may be different. Ours, using the primers reported by Gelernter et al (1999) are 376 and 419 bp. Heils et al (1996 ) report 484 and 528 bp; Wendland et al (2006) report 469 and 512 bp; and Nakamura et al (2000) refer to them as 14- and 16-repeat alleles. These are all the same.

## rs25531 SNP (MspI, L<sub>A</sub> / L<sub>G</sub>)

Hu et al (2006) reported that a SNP (rs25531, A/G) in the Long form of 5HTTLPR may have functional significance: The more common L<sub>A</sub> allele is associated with the reported higher basal activity, whereas the less common L<sub>G</sub> allele has transcriptional activity no greater than the S. These investigators suggest that in tests of association the L<sub>G</sub> alleles should be analyzed along with the S alleles.

The substitution of the G for A in the SNP, produces an MspI restriction site (CCGG) which forms the basis of the analysis strategy (Wendland et al, 2006). The sequence of the 5HTTLPR region was first reported by Heils, et al, and is reproduced below in the way they did to show the highly repetitive nature of the locus. There were two errors in their sequence in repeat units III and V (underlined) that have been corrected here. The forward and reverse primers we use now that yield amplicons of 376 and 419 bp, are shown in yellow highlight. The SNP, rs25531 is shown as "R" in green highlight and the second MspI site in repeat unit XIV is shown in teal highlight. The insertion/deletion is shown underlined in lower case in repeat units VI, VII, and VIII.


TCTCCCGCCTGGCGTTGCCGCTCTGAATGCCAGCACCTAACCCCTAATGT

I CCCTAC TGCA GCCCCCC AGCAT

II	CCCCC	TGCA	ACCTCCC	AGCA	
III	ACTCCC	TGTA	CCCCTCCT	AGGAT	
IV	CGCTCC	TGCA	TCCCCC	ATTATC	
V	CCCCC	TTCA	CCCCTCGC	GGCAT	
VI	CCCCC	TGCA	<u>cccc</u>	<u>R</u> gcat	R = A or G
VII	<u>cccccc</u>	<u>tgca</u>	<u>gccccccc</u>	<u>agcat</u>	
VIII	<u>ctcccc</u>	<u>tgca</u>	CCCCC	AGCAT	
IX	CCCCC	TGCA	GCCCTTCC	AGCA	
X	TCCCCC	TGCA	CCTCTCCC	AGGAT	
XI	CTCCCC	TGCA	ACCCCC	ATTAT	
XII	CCCCC	TGCA	CCCCTCGC	AGTAT	
XIII	CCCCC	TGCA	CCCCC	AGCATC	
XIV	CCCCCA	TGCA	CCC <u>CC</u>	<u>GG</u> CAT	
XV	CCCCC	TGCA	CCCCTCC	AGCAT	
XVI	TCTCCT	TGCA	CCCTACC	AGTAT	

TCCCCCGCATCCCGGCCTCCAAGCCTCCCGCCCACCTTGCGGTCCCCGCC

The forward primer has a fluorescent label attached to its 5' end. To analyze this SNP, the full length amplicons (from above) are incubated with the restriction enzyme MspI. The G allele which has the MspI restriction site (CCGG) will yield a product of 152 bp, whereas the A allele, which lacks the restriction site does not. A second MspI site 93 bp from the 3' end of the amplicon provides a positive control for the restriction reaction yielding cut products of 326 or 283 bp for the L and S alleles, respectively.



**5HTTLPR MspI restriction digest. The figure above is reproduced from a run from an ABI PRISM® 3130xl Genetic Analyzer. The amplicons are labeled with size in base pairs. The figure above shows all possible restriction fragments, and their source. The red peaks are size standards (Genescan ROX 500).**

To summarize, L<sub>G</sub> alleles yield 152 bp fragments, L<sub>A</sub> alleles yield 326 bp fragments and S alleles yield 283 bp fragments when incubated with MspI. For the following genotypes, the results for the PCR reaction, followed by MspI digest would be:

	<b>PCR</b>	<b>MspI</b>
L <sub>A</sub> /L <sub>A</sub>	419/419	326/326
L <sub>G</sub> /L <sub>G</sub>	419/419	152/152
L <sub>G</sub> /L <sub>A</sub>	419/419	152/326
S/S	376/376	283/283
S/L <sub>A</sub>	376/419	283/326
S/L <sub>G</sub>	376/419	283/152

You may notice that there is a third MspI site 30 bp from the end of the amplicon. This never shows up. Since the enzyme cleaves all of the sites equally well, only the smallest fragment with the 5' fluorescent label is visualized. If you were to run these cut products on an agarose gel, all could be visualized with ethidium bromide or other dye.

To analyze this SNP, the PCR products from above are used. After determining the genotype of the samples from above (e.g., LL, LS or SS), the PCR plate is prepared for MspI (#R106L, NEB, Ipswich, MA) restriction digest

### **Sample Preparation for MspI Digest**

95 °C 10 min  
65 °C 30 min  
4 °C hold

## MspI Restriction Digest Master Mix for 10 $\mu$ L reactions

(8  $\mu$ L Master mix + 2  $\mu$ L PCR product)

Component	1 Tube vol ( $\mu$ L)	100 Tubes vol ( $\mu$ L)	Concentration of component in:		
			Stock	Master Mix	Reaction
Water	6.8	680			
NEB buffer 2	1.0	100	10 x	0.125 x	1 x
MspI	0.2	20	20 units/ $\mu$ L	4 units	4 units
Total volume ( $\mu$ L)	8.0	800			

### MspI Digest Protocol

8  $\mu$ L of master mix + 2  $\mu$ L 5HTTLPR PCR product

37 °C 3 hours

65 °C 20 min

4 °C hold

### MspI Digest Electrophoresis

Size standard mixture for 100 samples:

500  $\mu$ L water

500  $\mu$ L Hi-Di formamide

25  $\mu$ L Genescan 500 Rox

Add 1  $\mu$ L digest product to 9  $\mu$ L size standard mix

Samples are analyzed on an ABI PRISM® 3130xl Genetic Analyzer using standard company protocols without modification

**Citation: When reporting results for this locus, please cite *Anchordoquy et al, 2003* as the analytical method used for genotyping.**

#### References:

Anchordoquy, H. C., McGeary, C., Liu, L., Krauter, K.S. and Smolen, A. (2003). Genotyping of three candidate genes following whole genome preamplification of DNA collected from buccal cells. *Behavior Genetics*. **33**: 73-78.

Don, R.H., Cox, P.T., Wainwright, B.J., Baker, K. and Mattick, J.S. (1992). "Touchdown" PCR to circumvent spurious priming during gene amplification. *Nucleic Acids Research* **19**: 4008.

Gelernter, J., Cubells, J.F., Kidd, J.R., Pakstis, A.J. and Kidd, K.K. (1999). Population Studies of Polymorphisms of the serotonin transporter protein gene. *American Journal of Medical Genetics (Neuropsychiatric Genetics)* **88**: 61–66.

Greenberg, B.D., Tolliver, T.J., Huang, S.J. Li, Q., Bengel, D., & Murphy D.L. (1999). Genetic variation in the serotonin transporter promoter region affects serotonin uptake in human blood platelets. *American Journal of Medical Genetics* **88**: 83-87.

Heils, A., A. Teufel, S. Petri, G. Stober, P. Riederer, D. Bengel, and K. P. Lesch. (1996). Allelic Variation of the Human Serotonin Transporter Gene Expression. *Journal of Neurochemistry* **66**: 2621-2624.

Hu, X., Oroszi, G., Chun, J., Smith, T.L., Goldman, D., and Schuckit, M.A. (2005). An expanded evaluation of the relationship of four alleles to the level of response to alcohol and the alcoholism risk *Alcoholism: Clinical and Experimental Research* **29**: 8-16.

Kaiser, R., Muller-Oerlinghausen, B., Filler, D., Tremblay, P. B., Berghofer, A., Roots, I., & Brockmoller, J. (2002) Correlation between serotonin uptake in human blood platelets with the 44-bp polymorphism and the 17-bp variable number of tandem repeat of the serotonin transporter. *American Journal of Medical Genetics* **114**: 323-328.

Lesch, K. P., Bengel, D., Heils, A., Sabol, S. Z., Greenberg, B. D., Petri, S., Benjamin, J., Muller, C. R., Hamer, D. H., and Murphy, D. L. (1996). Association of anxiety-related traits with a polymorphism in the serotonin transporter gene regulatory region. *Science* **274**:1527–1531.

Lichter, J.B., Barr, C.L., Kenedy, J.L., Van Tol, H.H. M., Kidd, K.K., & Livak, K.J. (1993). A hypervariable segment in the human dopamine receptor D4 (DRD4). *Human Molecular Genetics*, **2**: 767-773.

Nakamura, M., Ueno, S., Sano, A. and Tanabe, H. (2000). The human serotonin transporter gene linked polymorphism (5-HTTLPR) shows ten novel allelic variants. *Molecular Psychiatry* **5**: 32–38.

Nobile, M., Begni, B., Giorda, R., Frigerio, A., Marino, C., Molteni, M., Ferrarese, C., & Battaglia, M.J. (1999). Effects of serotonin transporter promoter genotype on platelet serotonin transporter functionality in depressed children and adolescents. *Journal of the American Academy of Child and Adolescent Psychiatry* **38**: 1396-1402.

Ramamoorthy, S., Bauman, A.L., Moore, K.R., Han, H., Yang-Feng, T., Chang, A.S., Ganapathy, V. and Blakely, R. D. (1993). Antidepressant- and cocaine-sensitive human serotonin transporter: molecular cloning, expression, and chromosomal localization. *Proceedings of the National Academy of Sciences USA* **90**: 2542–2546.

Wendland, J.R., Martin, B.J., Kruse, M.R, Lesch, K.-P and Murphy, D.L. (2006). Simultaneous genotyping of four functional loci of human SLC6A4, with a reappraisal of 5-HTTLPR and rs25531.

*Molecular Psychiatry* **11**: 224-226.

Willeit, M., Stastny, J., Pirker, W., Praschak-Rieder, N., Neumeister, A., Asenbaum, S., Tauscher, J., Fuchs, K., Sieghart, W., Hornik, K., Aschauer, H.N., Brucke, T. and Kasper, S. (2001) No evidence for in vivo regulation of midbrain serotonin transporter availability by serotonin transporter promoter gene polymorphism. *Biological Psychiatry* **50**: 8-12.

Williams, R.B., Marchuk, D.A., Gadde, K.M., Barefoot, J.C., Grichnik, K., Helms, M.J., Kuhn, C.M., Lewis, J.G., Schanberg, S.M., Stafford-Smith, M., Suarez, E.C., Clary, G.L., Svenson, I.K. and Siegler, I.C. (2003). Serotonin-related gene polymorphisms and central nervous system serotonin function. *Neuropsychopharmacology* **28**: 533-541.